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# Instructions for use Product group Ophthalmological Instruments

## Reusable surgically invasive instruments

No.	product group	No.	product group
1	01 Ophthalmological forceps	13	13 elevators
2	02 Ophthalmological knives	14	14 Ophthalmological Hooks
3	03 Marking instruments	15	15 chisels
4	04 calipers	16	16 loop instrument
5	05 General surgical scissors	17	17 eye spatulas / eye spoons
6	06 Ophthalmic scissors	18	18 eye probes
7	07 Callus scissors	19	19 eye retractors
8	08 Iris Scissors	20	20 lacrimal sac retractors
9	09 Enucleation scissors	21	21 needle holders
10	10 tenotomy scissors	22	22 bone punch
11	11 corneal trephine	23	23 Vascular clamp
12	12 terminal	24	24 scalpel handle

## Invasive reusable instruments

No.	product group	No.	product group
25	25 eye retractor	27	27 therapeutic contact lenses
26	26 tear duct probes	28	28 Eye Magnet

## Non-invasive reusable instruments

No.	product group	No.	product group
29	29 bandage scissors	--	--

Valid from: 2025-01-13

This instruction manual  
replaces the previous version:

1. Rev 01 / 2024-08-30

	<b>IMPORTANT PRODUCT INFORMATION PLEASE READ CAREFULLY BEFORE ANY CLINICAL USE!</b>
	The products are non-sterile and must undergo a processing process before clinical use.
	MEDICAL DEVICE IN THE SENSE OF REGULATION (EU) 2017/745 MDR

## Contents

<b>1</b>	<b>METAINFORMATICS</b>	<b>4</b>
1.1	MANUFACTURER	4
1.2	RELATED DOCUMENTS	4
1.3	SCOPE	4
1.3.1	Reusable surgically invasive instruments	4
1.3.2	Reusable invasive instruments	4
1.3.3	Reusable non-invasive instruments	4
<b>2</b>	<b>WHOLESALE MATERIALS</b>	<b>5</b>
<b>3</b>	<b>PRODUCT DESCRIPTION</b>	<b>5</b>
3.1	APPLICATION AREA	5
3.2	INTENDED PURPOSE & INDICATION	5
3.2.1	Reusable surgically invasive instruments	5
3.2.2	Reusable invasive instruments	7
3.2.3	Reusable non-invasive instruments	7
<b>4</b>	<b>CONTRAINDICATION</b>	<b>7</b>
4.1	GENERAL CONTRAINDICATIONS	7
4.2	PRODUCT-RELATED CONTRAINDICATIONS	7
<b>5</b>	<b>UNWANTED SIDE EFFECTS AND COMPLICATIONS</b>	<b>8</b>
5.1	GENERALLY	8
5.2	TREATMENT-RELATED	8
5.3	PRODUCT-RELATED	8
<b>6</b>	<b>WARNINGS</b>	<b>9</b>
<b>7</b>	<b>QUALIFICATION OF THE USER</b>	<b>10</b>
<b>8</b>	<b>INTENDED USER</b>	<b>10</b>
<b>9</b>	<b>CLINICAL BENEFITS</b>	<b>10</b>
<b>10</b>	<b>APPLICATION ENVIRONMENT</b>	<b>10</b>
<b>11</b>	<b>COMBINATION PRODUCTS &amp; ACCESSORIES</b>	<b>10</b>
11.1	SCALPELS WITH INTERCHANGEABLE BLADES	10
<b>12</b>	<b>APPLICATION &amp; SAFETY INSTRUCTIONS</b>	<b>10</b>
<b>13</b>	<b>STORAGE INSTRUCTIONS</b>	<b>11</b>
<b>14</b>	<b>PREPARATION</b>	<b>11</b>
14.1	WARNINGS	11
14.2	PLACE OF USE	11
14.3	TRANSPORT	11
14.4	PREPARATION FOR DECONTAMINATION	11
14.5	PRE-CLEANING	11
14.6	CLEANING / DISINFECTION	11
14.7	STERILIZATION	12
14.8	ADDITIONAL INFORMATION	12
<b>15</b>	<b>FUNCTIONAL TEST</b>	<b>12</b>
<b>16</b>	<b>MAINTENANCE AND CARE</b>	<b>13</b>

17	LIFETIME OF THE PRODUCTS.....	13
18	SERIOUS INCIDENTS .....	13
19	DISPOSAL.....	13
20	LIABILITY AND WARRANTY .....	13
21	SERVICE AND REPAIR.....	14
21.1	RETURN TRANSPORT.....	14
21.2	APPLICABLE STANDARDS – REFERENCES .....	14
22	SYMBOL EXPLANATIONS.....	15

# 1 METAinformatics

## Dear Customer!

By purchasing this instrument, you are receiving a high-quality product, the proper handling and use of which is described below. In order to keep risks to patients and users as low as possible, we ask you to read the instructions for use carefully and to follow them.

### 1.1 MANUFACTURER

**Carl Teufel GmbH & Co. KG**  
Tuttlinger Straße 30  
D- 78576 Emmingen - Liptingen  
Tel: +49 (0) 7465-1022  
Fax: +49 (0) 7465-1024  
<http://www.carl-teufel.de>  
office@carl-teufel.de

SINGLE REGISTRATION NUMBER

**DE-MF-000013978**

### 1.2 RELATED DOCUMENTS

Additional information may be required to use the products described here. This additional information, e.g. instructions for proper processing and sterilization, can be found on our homepage <http://www.carl-teufel.de> find.

<http://www.carl-teufel.de>

### 1.3 SCOPE

The scope of this instruction manual refers to all ophthalmic instruments that can be assigned to the respective product family and are shown in our catalogues:

#### 1.3.1 Reusable surgically invasive instruments

The conformity assessment procedure was carried out together with our notified body with the identification number 0483 with a focus on reprocessing. The products listed in the table are therefore marked with CE and the identification number.

product groups	product groups	product groups
01 Ophthalmological forceps	11 corneal trephine	21 needle holders
02 Ophthalmological knives	12 terminal	22 bone punch
03 Marking instruments	13 elevators	23 Vascular clamp
04 calipers	14 Ophthalmological Hooks	24 scalpel handle
05 General surgical scissors	15 chisels	--
06 Ophthalmic scissors	16 loop instrument	--
07 Callus scissors	17 eye spatulas / eye spoons	--
08 Iris Scissors	18 eye probes	--
09 Enucleation scissors	19 eye retractors	--
10 tenotomy scissors	20 lacrimal sac dilators	--

#### 1.3.2 Reusable invasive instruments

product groups	product groups	product groups
25 eye retractor	27 therapeutic contact lenses	--
26 tear duct probes	28 Eye Magnet	--

#### 1.3.3 Reusable non-invasive instruments

product groups	product groups	product groups
29 bandage scissors	--	--

## 2 WHOLESALE MATERIALS

- stainless steel
- titanium alloys
- Therapy contact lenses are made of silicone

## 3 PRODUCT DESCRIPTION

### 3.1 APPLICATION AREA

Our instruments are used in ophthalmology.

### 3.2 INTENDED PURPOSE & INDICATION

The instruments are intended for ophthalmological, surgical use with the intended purposes listed.

#### 3.2.1 Reusable surgically invasive instruments

##### 3.2.1.1 [01] (L1702 Tweezers)

Ophthalmological-surgical instrument for grasping and holding tissue and other material.

##### 3.2.1.2 [02] (Q020101 Ophthalmological knives )

Ophthalmic surgical instrument for making precise incisions in eye tissue during surgical procedures on the eye and adjacent structures.

##### 3.2.1.3 [03] (Q020203 Marking instruments)

Manual instrument for marking the cornea during ophthalmological surgery or preoperative treatment.

##### 3.2.1.4 [04] (Q0399 Caliper)

Hand-held ophthalmological measuring instrument consisting of two articulated legs for measuring diameter, length, angles and thicknesses of structures of the eye.

##### 3.2.1.5 [05] (L010499 Surgical scissors)

Scissors are hand-held, manual surgical instruments used for cutting body tissue, sutures or other materials.

##### 3.2.1.6 [06] (L010405 Ophthalmic scissors)

Hand-held ophthalmological-surgical instrument for cutting tissue, suture material during an ophthalmological procedure

##### 3.2.1.7 [07] (L01040501 Callus scissors)

Hand-held ophthalmic surgical instrument for cutting corneal tissue during surgical cataract extraction.

##### 3.2.1.8 [08] (L01040503 Iris scissors)

Hand-held ophthalmological-surgical instrument for cutting tissue during an operation on the iris of the eye (iris scissors), for cutting corneal tissue during surgical cataract extraction, for dissecting and/or separating tendon tissue in the eye or for cutting tissue in the anterior or posterior section during an operation on the eye.

##### 3.2.1.9 [09] (L01040502 Enucleation scissors)

Hand-held ophthalmic surgical instrument for cutting tissue during enucleation procedures on the eye and/or adjacent structures.

##### 3.2.1.10 [10] (L01040506 Tenotomy scissors and strabismus scissors )

Hand-held ophthalmological instrument for dissection and/or separation of tendon tissue in the eye in one procedure, usually to correct a deformity.

3.2.1.11 [11] (L170401 Corneal trephine )  
Ophthalmological-surgical knife for excising and removing round pieces of corneal tissue (corneal discs).

3.2.1.12 [12] (L170299 terminals)  
Surgical instrument for the temporary, atraumatic compression of an artery for the purpose of hemostasis (stopping or stopping bleeding) during an operation. \*

(\*excluded vessels: pulmonary arteries , aorta ascendens , arcus aortae , aorta descendens up to the bifurcation aortae , arteriae coronary arteries , arteries carotid communis , arteria carotid externa , arteria carotid interna , cerebral arteries , truncus brachiocephalicus , venae cordis , venae pulmonales, vena superior cava and vena inferior vena cava )

3.2.1.13 [13] (L0902 Elevator )  
Hand-held, manually operated surgical instrument for lifting/separating/stretching/widening (elevation) in ophthalmic surgery.

3.2.1.14 [14] (L1711 Ophthalmic hook)  
Hand-held ophthalmic surgical instrument for manipulating tissue structures in the eye or for removing foreign bodies from the eye.

3.2.1.15 [15] (L140602 chisel)  
hand-held surgical instrument for removing parts of bone or other hard tissue during orthopedic or plastic surgery.

3.2.1.16 [16] (L149008 Loop instruments)  
Hand-held ophthalmological surgical instrument for removing, among other things, the lens during enucleation.

3.2.1.17 [17] (L1706 Eye spatula/eye spoon)  
Handheld ophthalmic instrument for removing foreign bodies/objects that have penetrated or adhered to the surface of the eyeball.

3.2.1.18 [18] (L1705 Eye probes)  
A hand-held, rod-shaped, thin surgical instrument made of flexible metal with a blunt tip for examining or manipulating eye tissue or for dilating the tear duct.

3.2.1.19 [19] (L1701 Eye Wound Hook)  
Ophthalmic surgical instrument used to spread the edges of an incision in the eye and its adjacent structures during a surgical procedure.

3.2.1.20 [20] (L1701 Lacrimal Sac Spreader )  
Self-locking ophthalmic surgical instrument for exposure of the lacrimal sac

3.2.1.21 [21] (L0205 needle holder)  
Hand-held manual surgical instrument in microsurgery for grasping a sewing needle when it is pushed/pulled through tissue with the attached suture material during suturing.

3.2.1.22 [22] (L091302 Bone punches)  
Manual surgical instrument for taking soft tissue or bone samples during surgical procedures.

3.2.1.23 [23] (L0708010101 vascular clamp)  
Instrument for grasping, bringing together, compressing or holding an organ, vessel or tissue.\*

(\*excluded vessels: pulmonary arteries , aorta ascendens , arcus aortae , aorta descendens up to the bifurcation aortae , arteriae coronary arteries , arteries carotid communis , arteria carotid externa , arteria carotid interna , cerebral arteries , truncus brachiocephalicus , venae cordis , venae pulmonales, vena superior cava and vena inferior vena cava )

3.2.1.24 [24] (L010103 Scalpel handles )

Hand-held manual surgical instrument for interchangeable scalpel blades, with which the surgeon cuts tissue by hand.

### 3.2.2 Reusable invasive instruments

3.2.2.1 eye retractor

Surgical instrument used to retract the eyelids or other ocular tissues during ophthalmological examination, treatment, or surgical procedure.

3.2.2.2 Lacrimal duct probes:

A hand-held, rod-shaped, thin surgical instrument made of flexible metal with a blunt tip for examining or manipulating eye tissue or for dilating the tear duct.

3.2.2.3 therapeutic contact lenses:

A soft ophthalmic lens, usually made of a synthetic polymer/silicone, designed to be temporarily placed directly on the front surface of the eye (eyeball or cornea) for treatment and protection during surgical procedures on the eye. Protects the eye from damage during treatment.

3.2.2.4 Eye magnet:

Ophthalmological instrument with which metallic foreign bodies (usually containing iron or steel) in the eye tissue can be located and extracted.

### 3.2.3 Reusable non-invasive instruments

3.2.3.1 bandage scissors

A hand-held manual instrument for cutting bandages and fabric.

## 4 CONTRAINDICATION

### 4.1 GENERAL CONTRAINDICATIONS

- If the technical requirements for microsurgical intervention are not met immediately preoperatively or intraoperatively, the products listed here should not be used. By technical requirements we mean poor or inadequate lighting conditions.
- The instruments must not be used after contact with patients suffering from Creutzfeldt-Jakob disease (CJD or vCJD ) should not be reused or reprocessed.
- The instruments are contraindicated for any type of use outside of their intended use (see intended purpose).
- Application to the central cardiovascular system is contraindicated
- Due to their complexity, the instruments may only be used by trained and qualified personnel.

### 4.2 PRODUCT-RELATED CONTRAINDICATIONS

- None known

## 5 UNWANTED SIDE EFFECTS AND COMPLICATIONS

### 5.1 GENERALLY

- When using the products, unintentional damage to the surrounding tissue structures may occur, such as soft tissue injuries and tissue weakening and necrotization , even leading to fracture of bony structures.
- If the instruments are not prepared properly, inflammatory reactions and infections can occur.
- Furthermore, when handling instruments, injuries to the user can occur due to sharp-edged tips or pinching of tissue.

### 5.2 TREATMENT-RELATED

- perforation of the skin
- postoperative bleeding
- scars resulting from the surgical technique
- nodules on the eyelid edges
- Ectropion during eyelid surgery
- recurrence
- Skin shortage during eyelid surgery
- infections / inflammations
- Permanent visual impairment
- Rarely complete loss of vision
- seizures
- loss of consciousness

#### Special Complications of Blepharoplasty

- ptosis
- undercorrection
- overcorrection
- injury to the trochlea
- damage to the lacrimal gland
- blindness at 0.04%
- bleeding
- sunken lower eyelid
- drooping of the outer eyelid corner
- skin scars

#### Special Complications in Cataract Surgery

- an infection of the inside of the eye
- a swelling of the center of the retina
- a prolapse of the vitreous body due to injury to the posterior lens capsule
- secondary cataract

### 5.3 PRODUCT-RELATED

As part of our obligation to monitor the market, we have evaluated undesirable side effects and risks using comparable products from competitors that could also pose a risk to the respective product groups.

No.	product family	risks
1.	tweezers	<ul style="list-style-type: none"> <li>- Bruch</li> <li>- remaining pieces</li> <li>- injury to the surrounding area (tissue)</li> </ul>
2.	micro knife	<ul style="list-style-type: none"> <li>- metal particles in the eye</li> <li>- risk of blade breakage</li> <li>- cuts</li> </ul>
3.	corneal marker	None known
4.	calipers	None known
5.	surgical scissors	<ul style="list-style-type: none"> <li>- Bruch</li> </ul>

6.	Ophthalmic scissors	<ul style="list-style-type: none"> <li>- Possible unwanted perforation</li> <li>- remaining pieces</li> <li>- injury to the surrounding area (tissue)</li> <li>- injury to the user</li> </ul>
7.	corneal scissors	
8.	Iris Schere	
9.	enucleation scissors	
10.	STRABISMUS CORRECTION SCISSORS	
11.	tenotomy scissors	
12.	corneal trephine	None known
13.	artery clamp	<ul style="list-style-type: none"> <li>- Bruch</li> <li>- remaining pieces</li> <li>- injury to the surrounding area (tissue)</li> </ul>
14.	elevator	<ul style="list-style-type: none"> <li>- Bruch</li> </ul>
15.	hook, little hook	<ul style="list-style-type: none"> <li>- Bruch</li> <li>- deformation of the working ends</li> <li>- remaining pieces</li> <li>- injury to the surrounding area (tissue)</li> <li>- injury to the user</li> <li>- surgical extensions</li> </ul>
16.	chisel	None known
17.	snare instruments	None known
18.	spatula	None known
19.	eye probes	<ul style="list-style-type: none"> <li>- Bruch</li> <li>- deformation</li> <li>- loosening of components</li> </ul>
20.	lacrimal sac dilator	None known
21.	needle holder	<ul style="list-style-type: none"> <li>- loosening of components</li> <li>- remaining pieces</li> <li>- tissue injury</li> <li>- injury to the user</li> </ul>
22.	bone punches	None known
23.	vascular clamps	<ul style="list-style-type: none"> <li>- Bruch</li> <li>- remaining pieces</li> <li>- injury to the surrounding area (tissue)</li> </ul>
24.	scalpel handles	<ul style="list-style-type: none"> <li>- contamination residues from processing</li> </ul>
25.	eye retractor	None known
26.	lacrimal duct probes	None known
27.	therapeutic contact lenses	None known
28.	eye magnet	None known
29.	bandage scissors	None known

## 6 WARNINGS

### Danger!

- Use only for the specified purpose.
- All specifications described in the accompanying documentation must be strictly observed.
- Use only by qualified personnel
- Be careful when handling and transporting

### Danger!

The instruments of[Firma] are designed for surgical use only and must not be used for any other purpose. Improper handling and care as well as use for purposes other than those intended can lead to premature wear of the products.

### material incompatibility

The medical devices should not be used under any circumstances if the user or the personnel becomes aware that the patient has material intolerances.

**functional impairment**

Surgical instruments corrode and their function is impaired when they come into contact with aggressive substances. For this reason, it is essential to follow the processing and sterilization instructions.

**operating conditions**

To ensure the safe operation of the aforementioned products, correct maintenance and care is essential. In addition, a functional or visual inspection should be carried out before each use. For this reason, we refer you to the relevant sections in these instructions for use.

**storage**

There are no specific requirements for storing the products. However, we recommend storing the medical devices in a clean and dry environment.

## 7 QUALIFICATION OF THE USER

These products are only designed for qualified doctors or specialists with sufficient experience in the respective fields of application, e.g. ophthalmology. Explanations of detailed surgical procedures are therefore omitted.

## 8 INTENDED USER

These products are only intended for use by qualified physicians or specialists with sufficient experience in the respective fields of application, e.g. ophthalmology.

## 9 CLINICAL BENEFITS

The instruments are considered to be basic instruments in ophthalmology and the clinical benefit cannot be derived from the individual instrument. The clinical benefit is the optimization of treatment during a surgically invasive procedure. With the help of the instruments, surgical procedures can be carried out in a targeted manner, in accordance with the state of the art and recognized practice.

## 10 APPLICATION ENVIRONMENT

The products are intended for use in a controlled environment. The products must not be exposed to adverse/harsh environmental conditions. When used surgically, a sterile environment must be ensured.

## 11 COMBINATION PRODUCTS & ACCESSORIES

The products are not used with other products and are offered without accessories.

### 11.1 SCALPELS WITH INTERCHANGEABLE BLADES

Scalpels can be combined with blades according to DIN EN 27740. The scalpels are designed to be compatible with the figures according to DIN 58849-2.

## 12 APPLICATION & SAFETY INSTRUCTIONS

**delivery condition**

The medical devices are delivered in a non-sterile condition and must be prepared and sterilized by the user before the first and each subsequent use in accordance with the following instructions

## 13 STORAGE INSTRUCTIONS

Before use, products must be stored in an environment that preserves their packaging and purity. Dry atmosphere, no extreme temperatures, no exposure to sunlight, ionized radiation and contaminated particles. To avoid corrosion, pay special attention to the absence of chemicals in the immediate vicinity.

## 14 PREPARATION

### 14.1 WARNINGS

- Frequent reprocessing impairs the quality of the products.
- City water to be used must comply with COUNCIL DIRECTIVE 98/83/EC of 3 November 1998 on the quality of water intended for human consumption.
- This reprocessing instruction specifies the cleaning and disinfection agents used for validation. If an alternative cleaning agent and disinfectant (RKI or VAH listed) is used, the responsibility lies with the reprocessor.
- Reassemble disassembled products before sterilization.

### 14.2 PLACE OF USE

- The first steps of proper preparation begin in the operating room.
- Coarse soiling, residues of hemostatic agents, skin disinfectants and lubricants, as well as corrosive medicines, should be removed before the instruments are put down, if possible. Dry disposal (moistened, closed system) is preferred wherever possible. Avoid allowing residues to dry!
- Long waiting times until processing, e.g. overnight or over the weekend, should be avoided for both types of disposal (<6 hours).

### 14.3 TRANSPORT

- The products must be disposed of dry immediately after use. This means that the products must be transported wet in a closed container from the place of use to the processing site so that the products do not dry out.

### 14.4 PREPARATION FOR DECONTAMINATION

The products should, if possible, be disassembled before the subsequent processing steps or opened before being sent to the next processing steps. Rinsing shadows should be avoided. The products must be processed in suitable sieve baskets or rinsing bowls (select size according to product). The products should be fixed in the cleaning basket with a minimum distance between them. Overlapping should be avoided in order to rule out damage to the products during the cleaning process.

### 14.5 PRE-CLEANING

Products in an ultrasonic bath for 5 minutes using 0.5% Neodisher Pre-clean with Medizym. Rinse products under cold city water of drinking water quality (<40°C) until all visible dirt has been removed. Stubborn dirt should be removed with a soft brush or the cleaning mandrel provided. Cavities and lumina should be rinsed intensively (>30 sec) with cold city water of drinking water quality (<40°C) using a water pressure gun (or similar). Moving parts should be moved.

### 14.6 CLEANING / DISINFECTION

(washing machine, RDG):

- 1 minute pre-cleaning with cold city water Drinking water quality <40°C
- water drainage
- 3 minutes pre-cleaning with cold city water Drinking water quality <40°C
- water drainage
- 5 minutes cleaning at 55°C±5°C with 0.4% alkaline detergent ( Neodisher MediClean )
- water drainage
- 3 minutes neutralization (0.1% Neodisher ® Z) with cold city water Drinking water quality <40°C
- water drainage
- 2 minutes of rinsing with deionized water <40°C

The specific instructions of the cleaning machine manufacturer must be observed.

#### **Automatic disinfection**

Automatic thermal disinfection in cleaning and disinfection equipment, taking into account the national requirements for the A<sub>0</sub> value; e.g. A<sub>0</sub> value 3000:  
>5 minutes at 93°C±2°C with deionized water.

#### **Automatic drying**

Automatic drying according to the automatic drying process of the cleaning and disinfection device for at least 30 minutes (at 60°C±5°C in the washing chamber ). If necessary, subsequent manual drying with a lint-free cloth and blowing out lumens using sterile, oil-free compressed air.

## 14.7 STERILIZATION

Sterilization of the products using a fractionated pre-vacuum process (according to DIN EN ISO 17665-1) taking into account the respective national requirements. The products must be sterilized in suitable sterilization packaging.

Sterilization must be carried out using a fractionated pre-vacuum process with the following parameters:

- 134°C / 273.2°F,
- ≥5 minutes holding time,
- 3 pre-vacuum cycles

#### **Drying in vacuum for at least 20 minutes**

Flash sterilization is not suitable for products with lumens!

The autoclave manufacturer's instructions for use and the recommended guidelines for the maximum load of items to be sterilized must be observed. The autoclave must be installed, maintained, validated and calibrated correctly.

## 14.8 ADDITIONAL INFORMATION

The reprocessor is responsible for ensuring that the reprocessing actually carried out with the equipment, materials and personnel used in the reprocessing facility achieves the desired results. This usually requires validation and routine monitoring of the process and the equipment used.

## 15 FUNCTIONAL TEST

Inspect products after reprocessing and before sterilization for the following aspects:

- cleanliness
- Damage including, but not limited to, signs of corrosion (rust, pitting ), discoloration, deep scratches, peeling, cracks and wear.

- Proper functioning, including, among other things, sharpness of cutting tools, flexibility of flexible products, mobility of hinges/joints/box locks and moving parts such as handles and ratchets. Deformation of the products is not permitted.
  - Missing or removed (ground off) part numbers.
  - Do not use products that are not working properly, are defective, have excessive wear, or have defaced markings or part numbers that are missing or removed (ground off).
- Check products for flawless surfaces, correct assembly and functionality. Do not use severely damaged products, products with illegible markings, signs of corrosion or blunt cutting edges. Reassemble disassembled products before sterilization.

## 16 MAINTENANCE AND CARE

The products do not require maintenance. Medical white oil can be used to care for the instruments. Theraband can be used to check the functionality of scissors, see also table.

No.	article number	Description	Application
1.	100005347	Medical Fluid 20	For all joint instruments
2.	20030	Thera-Band exercise band	For all scissors (functional check)

## 17 LIFETIME OF THE PRODUCTS

The service life of the products depends essentially on the careful handling and preparation of the products by the user. Carl Teufel GmbH & Co. KG guarantees a service life of 2 years for reusable instruments when used carefully and prepared as intended.

If the user comes to the conclusion that the function is no longer available, the medical device must not be reprocessed and sterilized again. The END of the service life of the medical device has been reached.

The manufacturer validates a simulated repeated reprocessing for reusable surgically invasive instruments under worst case conditions. The reprocessing cycles are 100 repetitions and are intended to support microbiological safety. After completion of the validation, the results are included in detail in the instructions for use.

## 18 SERIOUS INCIDENTS

Any serious incident relating to the medical device must be reported to the manufacturer and to the competent authority in the Member State where the user and/or patient are located.

## 19 DISPOSAL

Disposal of the products, packaging materials and accessories must be carried out in accordance with applicable national regulations and laws.

## 20 LIABILITY AND WARRANTY

Carl Teufel GmbH & Co. KG, as the manufacturer, is not liable for consequential damages resulting from improper use or handling. This applies in particular to non-compliant use for the defined purpose or disregard of the processing and sterilization instructions. This also applies to repairs or modifications to the product that were carried out by unauthorized personnel of the manufacturer. These liability exclusions also apply to warranty services.

## 21 SERVICE AND REPAIR

Do not carry out any repairs or modifications to the product yourself. Only authorized personnel from the manufacturer are responsible for this and are intended to do so. If you have any complaints, claims or comments regarding our products, please contact us.

### 21.1 RETURN TRANSPORT

Defective or non-conforming products must have completed the entire remanufacturing process before being returned for repair/service.




### 21.2 APPLICABLE STANDARDS – REFERENCES

When cleaning, disinfecting and sterilizing, particular attention should be paid to the following sources:

1. AKI – Guideline “Instrument preparation done right”
2. RKI – Recommendation: “Requirements for hygiene during the processing of medical devices DIN EN 285 Large steam sterilizers
3. DIN EN 13060 Small steam sterilizers
4. DIN EN ISO 15883-1-3 Cleaning and disinfection equipment
5. DIN EN 868 / DIN EN ISO 11607-1 Packaging for medical devices to be sterilized in the terminal packaging
6. DIN EN 556-1 Sterilization of medical devices – Requirements for final packaging
7. DIN EN ISO 17664-1 Reprocessing of healthcare products - Information to be provided by the medical device manufacturer for the reprocessing of medical devices - Part 1
8. DIN EN ISO 17665 Sterilization methods – moist heat
9. DIN EN ISO 14937 Sterilization of health care products
10. DIN EN ISO 11737-1 Sterilization of medical devices – Microbiological process Part 1
11. DIN EN ISO 11737-2 Sterilization of medical devices – Microbiological process Part 2
12. DIN 58946-7 Sterilization, steam sterilizers
13. (AKI = Working Group on Instrument Reprocessing / RKI = Robert Koch Institute)
14. DIN EN ISO 96298-4 Medical instruments - Definitions, measurement methods and tests - Part 4: Functional control in the reprocessing cycle

## 22 SYMBOL EXPLANATIONS

	Manufacturer		Store in a dry place
	date of manufacture		article number
	unsterile		Protect from sunlight
	unique product identification	<b>SRN</b>	manufacturer's registration number in the EUDAMED database
	Danger		Unique Device Identification, code for identifying a product
	Follow the instructions for use		batch number
	Federal law restricts this device to sale by or on the order of a physician, dentist, or licensed professional.		order number
	CE marking with identification number of the notified body.		medical device
	indication of the unit of measure		CE mark

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